




Self-Retaining Retractors

- Accommodates up to .062 K-Wire.
- Can be used as a distractor and/or retractor.
- Ratchet graduated up to 35mm to accurately assess the width for bone graft requirements.
- Useful in many types of procedures, such as fusions, fractures and nonunions.
- Narrow tips allow for easy entrance into small incisions.
- Serrated tips for secure positioning.

DELTA EXCLUSIVE


Modified Inge Retractor with K-Wire Guides

D 544 101
6-1/2" 35mm Graduated Ratchet



Inge Retractor

D 544 070
6-1/2" with Teeth
D 544 080
6-1/2" without Teeth



Modified Inge Retractor

D 544 100
6-1/2" 35mm Graduated Ratchet