


Super-Cut Scissors



Super-Cut Iris Scissors

Sharp Points


- D 416 400 3-1/2" Straight, Delicate
- D 416 420 3-1/2" Curved, Delicate
- D 416 440 4-1/2" Straight
- D 416 500 4-1/2" Curved

Super-Cut Stevens Scissors

Blunt Points

- D 416 740 4-1/2" Straight
- D 416 800 4-1/2" Curved
- D 416 840 5-1/4" Curved


Fine sharp blades for cutting delicate tissue during plastic, hand, and minor vascular surgery.



Super-Cut Operating Scissors

Sharp/Blunt Points



- D 414 100 5-1/2" Straight
- D 414 400 5-1/2" Curved



Super-Cut Jamison Scissors

- D 416 860 6-1/4" Curved, Blunt Points
- Very Delicate blunt scissors with long handles for fine tissue dissection



Operating and Dissecting Scissors



Iris Scissors


Sharp Points

- D 592 540 3-1/2" Straight
- D 592 560 3-1/2" Curved
- D 592 580 4-1/2" Straight
- D 592 620 4-1/2" Curved



Operating Scissors



- D 395 100
5-1/2" Straight, Sharp/Blunt Points
- D 395 120
5-1/2" Straight, Sharp/Sharp Points



Mayo Scissors

Beveled Blades, Blunt Points

- D 404 460 5-1/2" Straight
- D 404 480 5-1/2" Curved
- D 404 500 6-3/4" Straight
- D 404 520 6-3/4" Curved



Metzenbaum Scissors

Blunt Points

- D 406 480 6" Straight
- D 406 520 6" Curved
- D 406 560 7" Straight
- D 406 600 7" Curved