





Rongeurs



Micro Friedman Rongeur

1.5x5mm Bite



- D 534 300 5-1/2" Straight
- D 534 320 5-1/2" 30° Angled
- D 534 340 5-1/2" 90° Angled
- D R11 550 Replacement Springs (2)



Friedman Rongeur

4x11mm Bite


- D 534 520
- 5-1/2" Curved Delicate



Mini-Friedman Rongeur

2x8mm Bite, Delicate

- D 534 530 6" Straight
- D 534 540 6" Curved






Dean (Blumenthal) Rongeur




4x11mm Bite

- D 534 560 6" 30° Curve



Rongeurs



Lempert Rongeur
3x7mm Bite, Delicate
D 534 480 6" Straight
D 534 500 6" Curved





Lempert Rongeur
3x7mm Bite, Delicate
D 534 440 7-1/2" Straight
D 534 460 7-1/2" Curved



Stille-Luer Rongeur
10x18mm Bite
D 536 480 8-1/2" Straight
D 536 500 8-1/2" Curved


Rongeurs and Trephines



Leksell Rongeur

9", 15° Curve

- | | | | |
|-----------|-------------|-----------|-------------|
| D 536 820 | 3x16mm bite | D 536 780 | 6x16mm bite |
| D 536 760 | 4x16mm bite | D 536 790 | 7x16mm bite |
| D 536 770 | 5x16mm bite | D 536 800 | 8x16mm bite |



*Designed for bone biopsy.
Also facilitates removal of
broken bone screws.*



Michele Trephine

- D 544 290
7" Small, 3mm inside diameter, 5mm outside diameter
- D 544 300
7" Medium, 4mm inside diameter, 6mm outside diameter
- D 544 310
7" Large, 6mm inside diameter, 8mm outside diameter