


# Steinmann Pins

## Steinmann Pin Implants

- Precision manufactured in Germany from certified implant stainless steel.
- Smoothly hand polished, causing less trauma when removed.
- Steinmann Pins are sold six (6) to the package.




### 9" Smooth, Single Trocar Round End

D 892 030	5/64" (2.0mm) diameter
D 892 430	3/32" (2.4mm) diameter
D 892 830	7/64" (2.8mm) diameter
D 893 230	1/8" (3.2mm) diameter
D 893 530	9/64" (3.5mm) diameter
D 894 030	5/32" (4.0mm) diameter
D 894 530	.177" (4.5mm) diameter



### 9" Smooth, Single Diamond Round End

D 892 010	5/64" (2.0mm) diameter
D 892 410	3/32" (2.4mm) diameter
D 892 810	7/64" (2.8mm) diameter
D 893 210	1/8" (3.2mm) diameter
D 893 510	9/64" (3.5mm) diameter
D 894 010	5/32" (4.0mm) diameter
D 894 510	.177" (4.5mm) diameter



### 9" Smooth, Double End Trocar

D 892 040	5/64" (2.0mm) diameter
D 892 440	3/32" (2.4mm) diameter
D 892 840	7/64" (2.8mm) diameter
D 893 240	1/8" (3.2mm) diameter
D 893 540	9/64" (3.5mm) diameter
D 894 040	5/32" (4.0mm) diameter
D 894 540	.177" (4.5mm) diameter



### 9" Smooth, Double End Diamond

D 892 020	5/64" (2.0mm) diameter
D 892 420	3/32" (2.4mm) diameter
D 892 820	7/64" (2.8mm) diameter
D 893 220	1/8" (3.2mm) diameter
D 893 520	9/64" (3.5mm) diameter
D 894 020	5/32" (4.0mm) diameter
D 894 520	.177" (4.5mm) diameter

*Please inquire about availability of any size and style of Pin Implants not shown in this catalog.*

*Drawings shown on this page are not actual size.*