







K-Wire Driver and Gratloch™ Wire Bender



- Drive, bend and pull K-Wire, all in one tool
- One handed operation
- Powerless, ideal for office or hospital use
- Contains depthscale
- No set-up or maintenance
- Fully autoclavable
- For wire from .028" to .062"

K-Wire Driver

D 527 300
5-1/2" German Stainless




The Gratloch™ Wire Bender stabilizes and bends K-Wire simultaneously. Wire can be bent close to the bone single-handedly without fracture or osteotomy stability.





P.I. Gratloch™ Wire Bender

D 545 600
7-1/2" Hardened, Wear-resistant Stainless Steel
Maximum capacity .062" (1.6mm)
U.S. Patent #5,314,431

Wire Pullers and Benders



Wire Pulling Forceps
 D 528 270
 6-1/2"
 Maximum capacity 1mm



Narrow Nose Wire Pulling Forceps
 D 544 430
 7-1/2"
 3mm tip



K-Wire Bender
 D 559 920
 5-1/2" 3 Prong
 Bends up to 1.6mm (.062") K-Wire to 90° angle

